

Omega Laboratories Oral Fluid Drug Panels & Cutoff Levels

	Screening (Dilute)	Confirmation (Neat)	10 Panel <i>S10P</i>	6 Panel <i>S6P</i>
Amphetamines	100 ng/ml		●	●
Amphetamine		50 ng/ml	●	●
Methamphetamine / MDMA	40 ng/ml		●	●
Methamphetamine		50 ng/ml	●	●
MDA		50 ng/ml	●	●
Cocaine / Metabolite	5 ng/ml		●	●
Cocaine		8 ng/ml	●	●
Benzoylcegonine		8 ng/ml	●	●
Marijuana (THC)	1 ng/ml		●	●
THC		2 ng/ml	●	●
Opiates	10 ng/ml		●	●
Codeine		40 ng/ml	●	●
Hydrocodone		40 ng/ml	●	●
Phencyclidine	1 ng/ml		●	●
Phencyclidine (PCP)		10 ng/ml	●	●
Oxycodone	10 ng/ml		●	
Oxycodone		40 ng/ml	●	
Methadone	5 ng/ml		●	
Methadone		15 ng/ml	●	
Benzodiazepines	1 ng/ml		●	
Alprazolam		5 ng/ml	●	
Diazepam		5 ng/ml	●	
Barbiturates	20 ng/ml		●	
Secobarbital		60 ng/ml	●	
Pentobarbital		60 ng/ml	●	
Amobarbital		60 ng/ml	●	
Separate Panels				
Cotinine (Nicotine Metabolite)				
Cotinine	10 ng/ml	10 ng/ml	(Offered separate from 6 and 10 panels)	

Neat vs. Dilute: Screening cutoff levels listed as Dilute to account for device buffer. Confirmation listed as Neat in accordance with SAMHSA guidelines. For Neat screening level, multiply Dilute number times 3.

Validity Testing: All oral fluid panels include IgG human validity testing