

Omega Laboratories Oral Fluid Drug Panels & Cutoff Levels

	Screening (Neat)	Confirmation (Neat)	10 Panel <i>S10P</i>	6 Panel <i>S6P</i>
Amphetamines	300 ng/ml		●	●
Amphetamine		50 ng/ml	●	●
Methamphetamine / MDMA	120 ng/ml		●	●
Methamphetamine		50 ng/ml	●	●
MDA		50 ng/ml	●	●
Cocaine / Metabolite	15 ng/ml		●	●
Cocaine		8 ng/ml	●	●
Benzoylcegonine		8 ng/ml	●	●
Marijuana (THC)	3 ng/ml		●	●
THC		2 ng/ml	●	●
Opiates	30 ng/ml		●	●
Codeine		40 ng/ml	●	●
Hydrocodone		40 ng/ml	●	●
6-MAM (Heroin)		4.0 ng/ml	●	●
Phencyclidine	3 ng/ml		●	●
Phencyclidine (PCP)		10 ng/ml	●	●
Oxycodone	30 ng/ml		●	
Oxycodone		40 ng/ml	●	
Methadone	15 ng/ml		●	
Methadone		15 ng/ml	●	
Benzodiazepines	3 ng/ml		●	
Alprazolam		5 ng/ml	●	
Diazepam		5 ng/ml	●	
Barbiturates	60 ng/ml		●	
Secobarbital		60 ng/ml	●	
Pentobarbital		60 ng/ml	●	
Amobarbital		60 ng/ml	●	
Separate Panels			(Offered separate from 6 and 10 panels)	
Cotinine (Nicotine Metabolite)	30 ng/ml			
Cotinine		10 ng/ml		

Neat vs. Dilute: Screening and Confirmation listed as Neat in accordance with recent SAMHSA guidelines.

Validity Testing: All oral fluid panels include IgG human validity testing