

Ω Omega Laboratories Oral Fluid Drug Panels & Cutoff Levels

Oral Fluid 6-Panel (S6P)	Screening Cutoff (Dilute)	Confirmation Cutoff (Neat)
Amphetamine <i>Amphetamine (Adderall, Benzedrine, Dexedrine, Vyvanse)</i>	100 ng/ml	50 ng/ml
Methamphetamines <i>Methamphetamine, MDM, MDA, MDEA</i>	40 ng/ml	50 ng/ml 50 ng/ml
Cocaine <i>Cocaine, Benzoylcegonine</i>	5 ng/ml	8 ng/ml 8 ng/ml
Opiates <i>Codeine, Morphine, Hydrocodone (Vicodin, Lortab/Lorcet), 6-MAM (Heroin)</i>	10 ng/ml	40 ng/ml 40 ng/ml 40 ng/ml 4 ng/ml
Phencyclidine <i>PCP</i>	1 ng/ml	10 ng/ml
Marijuana (THC) <i>THC</i>	1 ng/ml	2 ng/ml

Oral Fluid 10-Panel (S10P) (In Addition to 6-Panel)	(Dilute)	(Neat)
Oxycodone <i>Oxycodone (OxyContin, Percodan/Percocet)</i>	10 ng/ml	40 ng/ml
Methadone <i>Methadone</i>	5 ng/ml	15 ng/ml
Benzodiazepines <i>Diazepam, Nordiazepam (Valium), Alprazolam, Alpha-hydroxyalprazolam (Xanax)</i>	1 ng/ml	5 ng/ml 5 ng/ml
Barbiturates <i>Secobarbital (Seconal), Butalbital (Fioricet, Fiorinal), Pentobarbital (Nembutal), Phenobarbital (Luminal), Amobarbital (Amytal, Tuinal)</i>	20 ng/ml	60 ng/ml 60 ng/ml 60 ng/ml 60 ng/ml 60 ng/ml

Separate Panels (Offered separate from 6 and 10 panels)	(Dilute)	(Neat)
Cotinine (Nicotine Metabolite) (SCOT) <i>Cotinine</i>	10 ng/ml	10 ng/ml

Neat vs. Dilute: Screening cutoff levels listed as Dilute to account for device buffer. Confirmation listed as Neat in accordance with SAMHSA guidelines

Validity Testing: All oral fluid panels include IgG human validity testing