

# Omega Laboratories Oral Fluid Drug Panels & Cutoff Levels

	Screening (Dilute)	Confirmation (Neat)	10 Panel <i>S10P</i>	6 Panel <i>S6P</i>
<b>Amphetamines</b>	100 ng/ml		●	●
Amphetamine		50 ng/ml	●	●
<b>Methamphetamine /MDMA</b>	40 ng/ml		●	●
Methamphetamine		Ecstasy (MDMA) 50 ng/ml	●	●
MDA		MDEA 50 ng/ml	●	●
<b>Cocaine /Metabolite</b>	5 ng/ml		●	●
Cocaine		8 ng/ml	●	●
Benzoylcegonine		8 ng/ml	●	●
<b>Marijuana (THC)</b>	1 ng/ml		●	●
THC		2 ng/ml	●	●
<b>Opiates</b>	10 ng/ml		●	●
Codeine		Morphine 40 ng/ml	●	●
Hydrocodone		6-MAM (Heroin) 40 ng/ml	●	●
<b>Phencyclidine</b>	1 ng/ml		●	●
Phencyclidine (PCP)		10 ng/ml	●	●
<b>Oxycodone</b>	10 ng/ml		●	
Oxycodone		40 ng/ml	●	
<b>Methadone</b>	5 ng/ml		●	
Methadone		15 ng/ml	●	
<b>Benzodiazepines</b>	1 ng/ml		●	
Alprazolam		Nordiazepam 5 ng/ml	●	
Diazepam		5 ng/ml	●	
<b>Barbiturates</b>	20 ng/ml		●	
Secobarbital		Butalbital 60 ng/ml	●	
Pentobarbital		Phenobarbital 60 ng/ml	●	
Amobarbital		60 ng/ml	●	
<b>Separate Panels</b>				
<b>Cotinine (Nicotine Metabolite)</b>				
Cotinine	10 ng/ml	10 ng/ml	(Offered separate from 6 and 10 panels)	

**Neat vs. Dilute:** Screening cutoff levels listed as Dilute to account for device buffer. Confirmation listed as Neat in accordance with SAMHSA guidelines

**Validity Testing:** All oral fluid panels include IgG human validity testing