

The iScreen™ Instant Saliva Screen from OMEGA LABORATORIES, INC.

The iScreen instant saliva drug screen offers you ease of collection and laboratory confirmation if required. The iScreen is a fast, donor-friendly test that provides results in just 10 minutes.


Advantage of iScreen over urine testing

1. No urine handling
2. No “shy bladder” delays
3. Very difficult to adulterate
4. Test multiple donors at the same time
5. Sample substitution virtually impossible
6. Least invasive collection
7. Short test time
8. More immediate detection time




Here's how iScreen works

1. Insert sponge end of collector into mouth
2. Make sure sponge is softened (3 minutes)
3. Insert collector in reservoir and twist to lock
4. Wait 10 minutes and read results



Advantages over laboratory testing

1. Instant results
2. Lower costs
3. Hire on the spot
4. Instant return to duty
5. No collection site needed
6. No travel to collection site
7. Virtually impossible to substitute
8. Minimal shipping costs



Detects Up to 6 Drugs

Phencyclidine	Cocaine
Amphetamines	Opiates (Morphine)
Marijuana	Methamphetamine

Ω OMEGA
Laboratories

400 N. Cleveland Avenue
Mogadore, Ohio 44260
www.omegalabs.net

800.665.5569 330.628.5748