

The iScreen from OMEGA LABORATORIES, INC.

The iScreen drug screen is an on-site (rapid) urine test offering low cost and convenient collection. The iScreen is used by removing the end cap and inserting the detection device into the collection cup. Once the detection device is inserted, the tests will begin to run. The test results will be visible in a few minutes and if no drugs are found the corresponding panel will develop a red line below the Control Line. Optional confirmation kits are available.


✓ Instructions for reading test results

✓ Drug name abbreviation on test panel for easy identification

✓ Control lines for test

✓ Test lines appear here to indicate no drugs were found

✓ End cap to cover test strips



Detects Up to 9 Drugs

Phencyclidine	Methamphetamine
Amphetamines	Benzodiazepine
Marijuana	Barbiturates
Cocaine	Methadone
Opiates (Morphine)	

 **OMEGA**
Laboratories

400 N. Cleveland Avenue
Mogadore, Ohio 44260
www.omegalabs.net
800.665.5569 330.628.5748

Copyright 2002, Omega Laboratories, Inc.